

# Bathing Cushion

The Bathing Cushion is used to help someone who has difficulty getting in and out of the bath



"Warm water reduces musculoskeletal pain."  
Arthritis Foundation



## Technical Specifications

Max user weight	150 kg / 24 stone
Max height	40 cm / 15 ¾ "
Min height	Flat
Seat depth	44 cm / 17 ½ "
Width	58 cm / 23 "
Total Product weight	2 kg / 4 ½ lbs
Back/headrest length	33 cm / 13 ¾ "
Airflo Compressor	Mk3

Product code: MPCA030200

### Easy to fit

- Fitted and removed within seconds
- Some air can be left in the backrest for greater levels of comfort
- The design adapts to the size and shape of the bath, so no expensive bathroom alterations are required

### Safe to use

- Only lowers a person into the bath if there is enough battery power to raise them back out again
- Battery operated – no electric cables
- Fully waterproof hand control (IP67)

### Simple to maintain

- Lightweight and portable, it can be removed and cleaned with ease
- Rinse with clean water after each bath and leave it to dry inflated
- Compact and easy to store

### Improves wellbeing

- Perfect for a person who has difficulty getting in and out of the bath
- The cushion deflates allowing the user to lie flat in the bath, fully covered in water so benefiting from the therapeutic properties of bathing
- The user can take the bathing cushion with them when travelling away
- Can be removed from the bath when not in use so enabling a more flexible environment for shared bathrooms

